V. Mueller™

en	Neuro/Spine Products Spetzler Headrest System
fr	Produits neuro-rachidiens Système de repose-tête Spetzler
de	Neurologische/Wirbelsäulen-Produkte Spetzler-Kopfstützensystem
it	Prodotti neuro/spinali Sistema poggiatesta Spetzler
es	Productos neurológicos/medulares Sistema de reposacabezas Spetzler
pt	Produtos Neuro/Spine Sistema de Fixação Craniana Spetzler
nl	Neuro-/wervelkolomproducten Spetzler hoofdsteunsysteem
sv	Neuro/Spine Products Spetzler Headrest System
da	Neuro/Spine Produkter Spetzler Hovedstøtte System
fi	Neuro-/selkärätkatuotteet Spetzler-päätukijärjestelmä
no	Nevrologi/spinal-produkter Spetzler hodestøttesystem
el	Νευρολογικά/Σπονδυλικά προϊόντα Σύστημα στηρίγματος κεφαλής Spetzler
ja	Neuro/Spine 製品 Spetzler ヘッドレストシステム
cs	Neurologické / spinální produkty Systém opěrky hlavy Spetzler
pl	Produkty do zabiegów neurochirurgicznych i na kręgosłupie Wspornik głowy Spetzler
tr	Sinir / Omurga Ürünleri Spetzler Baş Dayanağı Düzeni
ro	Produse neuro-vertebrale Sistem de sprijin pentru cap Spetzler

USA Rx Only

Snowden-Pencer is a trademark and/or registered trademark of CareFusion Corporation, or one of its subsidiaries.

© 2012, CareFusion Corporation or one of its subsidiaries. All rights reserved.

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


Device Components

Figure 1: M-1525 Adjustable Base Unit

- (a) Transitional Member
- (b) Attachment Rods
- (c) End Bracket
- (d) Hex Wrench
- (e) Connecting Tube
- (f) Tension Adjustment Knob
- (g) Locking Lever
- (h) Socket Screw
- (i) Nylon Retaining Screw

Figure 2: M-1570 Swivel Adapter

- (a) Knob
- (b) Retaining Ring
- (c) Knob
- (d) Dowel Pin
- (e) Swivel Adapter Base

Figure 3

- (a) Skull Clamp Threaded Hole
- (b) Knob
- (c) Swivel Adapter
- (d) Threaded Screw

Figure 4: M-1500 Skull Clamp

- (a) Force Gauge
- (b) Torque Screw
- (c) Ratchet Arm
- (d) Locking Ratchet
- (e) Starbursts
- (f) Skull Clamp Base
- (g) Index Knob
- (h) Rocker Arm
- (i) Skull Pin

Figure 5

- (a) Skull Pin Receptacle
- (b) Threaded Area
- (c) Pressure Gauge Knob

Figure 6

- (a) Correct
- (b) Not Correct

Indications For Use

V.Mueller Neuro/Spine Products comprising surgical instruments that provide rigid 3-point skeletal fixation for the supine, prone, lateral or sitting positions, craniotomies in all positions and cervical laminectomies (upper dorsal approach).

How Supplied

The Spetzler Headrest System is packaged as non-sterile. Skull pins are available sterile or non-sterile.

Limitations on Reprocessing

Repeated reprocessing has minimal effect on these devices. End of life is normally determined by wear and damage due to use. Spetzler Disposable Skull Pins are for single use only and should be disposed after use.

Warnings

Devices shall be used in accordance with these instructions for use. Read all sections of this insert prior to use. Improper use of this device may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and cause a serious injury to the patient or health care provider.

Any variation between what is suggested here and your facility's policies and procedures should be brought to the attention of the appropriate responsible persons for direction before using your instruments. Prior to use, all instruments should be inspected to insure proper function and condition. Do not use instruments if they do not perform satisfactorily in an operational test.

Do not position skull pins over the frontal sinuses, the temporal fossa, and major vessels and nerves of the skull. Avoid puncture of the superficial temporal artery.

Special care should be exercised on pediatric patients due to thin skull consistency. Spetzler Skull Clamp should not be used on children under five (5) years of age.

Special care should be exercised with use on the elderly, on restored surgical areas surrounding previous burr holes, or any diseased bone because of varying consistencies and thickness from that of healthy bone.

A minimum of 60 lbs force must be applied to the skull clamp. Skull clamp force above 80 pounds could cause harm to the patient. **Do Not Exceed 80lbs of force when turning the Torque Knob.**

Skull pins should penetrate the scalp at an angle that intersects the center of the cranial vault. Never sharpen skull pins. Dull or damaged pins should be discarded and replaced.

It is not recommended, and may not be necessary, to sterilize the Spetzler Headrest System. If sterilization is indicated, Ethylene Oxide sterilization may be used. **Do not steam sterilize the headrest system.**

The Spetzler Skull Pin must be fully seated so that there are no gaps between the shoulder of the skull pin and the receptacle into which it has been inserted. Skull pins should not be centered over the frontal sinuses or the temporal area of the skull. Failure to follow these warnings may result in serious patient injury or death. To maintain proper fixation, skull pins should be of equal length.

Failure to properly position the patient and to fully tighten all adjustment portions of the device may result in patient injury.

Cautions

If there are any variations between these instructions and either your facility's policies and/or your cleaning/sterilization equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your devices.

Use of device for a task other than that for which it is intended will usually result in a damaged or broken device.

Prior to use, inspect device to ensure proper function and condition. Do not use devices if they do not satisfactorily perform their intended function or have physical damage.

Pay special attention to all threads; be sure there is no thread damage. Do not use if threads are damaged or not uniform in shape.

The Spetzler Headrest System is supplied non-sterile and must be properly draped prior to each use.

Visually inspect the drape to ensure that no holes and/or tears occurred during setup. If a barrier breach is found, discard drape and repeat the draping procedure.

Ensure that all parts of the device in the sterile field are properly draped or sterile.

Use only neutral pH (6-8) detergent solutions.

Pre-processing Instructions

Cleaning and sterilization is recommended in the event that the device is contaminated during use.

Initiate cleaning of device within 2 hours of use.

Transport devices via the institution's established transport procedure.

Remove excess gross soil as soon as possible after use by rinsing or wiping the device.

Turn the index knob to the unlocked position.

Disassembly is required prior to cleaning and sterilization.

All devices must be processed in the completely open position (i.e. brackets, levers, etc.) to allow solution contact of all surfaces.

Manual Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Prepare the enzymatic / neutral pH detergent solution, utilizing tap water with a temperature range of 27°C to 44°C (81°F to 111°F), per manufacturer's instructions.
3. Place device in the open/relaxed position, and completely immerse in the detergent solution and allow device to soak for a minimum of 5 minutes. Actuate all movable parts during the initiation of the soak time.
4. Using a soft bristled brush, remove all visible soil from the device. Actuate device while brushing, paying particular attention to hinges, crevices and other difficult to clean areas. **Note:** It is recommended that the detergent solution is changed when it becomes grossly contaminated (bloody and/or turbid).
5. Rinse the device by completely immersing in tap water with a temperature range of 27°C to 44°C (81°F to 111°F), for a minimum of 30 seconds to remove any residual detergent or debris.
6. Dry the device with a clean, lint-free towel.
7. Visually examine each device for cleanliness.
8. If visible soil remains, repeat cleaning procedure.

Automatic Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Clean the devices via the automatic cleaning parameters below.

PHASE	MINIMUM RECIRCULATION TIME	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION (IF APPLICABLE)
Pre-Wash 1	15 Seconds	Cold Tap Water 1°C - 16°C (33°F - 60°F)	N/A
Enzyme Wash	1 Minute	Hot Tap Water 43°C - 82°C (110°F - 179°F)	• Detergent: pH-neutral/enzymatic detergent • Concentration: Per the detergent manufacturer's recommendations
Wash 1	2 Minutes	Tap Water 43°C - 82°C (110°F - 179°F)	• Detergent: pH-neutral cleanser • Concentration: Per the detergent manufacturer's recommendations
Rinse 1	15 Seconds	Tap Water 43°C - 82°C (110°F - 179°F)	N/A
Pure Rinse	10 Seconds	Purified Water 43°C - 82°C (110°F - 179°F)	N/A
Drying	N/A	N/A	N/A

3. If visible moisture is present dry the instrument with a clean, lint-free towel.
4. Visually examine each instrument for cleanliness.
5. If visible soil remains, repeat cleaning procedure.

Inspection/Maintenance

Proper care and handling is essential for satisfactory performance of any surgical device. The previous cautions should be taken to ensure long and trouble-free service from all your surgical devices. Inspect devices before each use for broken, cracked, tarnished surfaces, movement of hinges, and chipped or worn parts. If any of these conditions appear, do not use the device. Return devices to an authorized repair service center for repair or replacement.

Check for smooth movement of knobs and rocker arm without excessive "play". Locking mechanisms should be free of nicks and present a continuous edge.

Before sterilizing, lubricate the device with instrument milk or a steam permeable/ water soluble lubricant, following the lubricant manufacturer's instructions.

Let devices drip dry for three (3) minutes before packaging for sterilization.

Packaging

Devices can be loaded into dedicated packaging systems. Sterilization wrap material must be cleared for the applicable sterilization modality by your country's regulatory body. Use in accordance with packaging manufacturer's sterilization instructions being sure to protect jaws and cutting edges from damage.

Sterilization

Sterilization of instruments may be accomplished by Ethylene Oxide.

Do Not Steam Sterilize the Headrest System.

Check to ensure index knob is in the unlocked position prior to sterilization.

All devices must be processed in the completely open position (i.e. brackets, levers, etc.) to allow sterilant contact of all surfaces.

All devices with concave surfaces shall be configured so that water pooling does not occur.

EO Sterilization Parameters

Sterilizer Cycle: 100% Ethylene Oxide (EO)

Minimum Preconditioning Time: 30 minutes

Minimum EO Gas Concentration: 725 mg/L

Minimum Temperature: 55°C (130°F)

Minimum Humidity: 50-80%

Minimum Exposure Time: 60 minutes

Minimum Aeration: 8 hours @ 43°C (110°F)

Sterilization Configuration: Wrapped (2 layer 1-ply or 1 layer 2-ply)

Storage

After sterilization, devices must remain in sterilization packaging and be stored in a clean, dry environment.

The Base Unit must be stored in the locked position with the Transitional Member inserted. Failure to lock both components will result in damage to the Transitional Member receptacle.

Assembly Instructions

Attach the Adjustable Base Unit by inserting the Attachment Rods into the operating room table. The Left End Bracket is adjustable and may be repositioned to adapt the Base Unit to varying operating room table configurations. To adjust End Bracket, remove Hex Wrench from its position by turning it counterclockwise. Insert Hex Wrench into the socket screw located at the base of the Left Adjustable End Bracket. Turn Hex Wrench clockwise to loosen End Bracket and slide to desired position. Tighten socket screw until End Bracket is securely in place. **Figure 1**

Insert threaded screw located on the smaller Swivel Adapter starburst into the threaded hole on the Adjustable Base Unit's Transitional Member. Turn knob clockwise to tighten. Insert threaded screw located on the larger Swivel Adapter starburst to the desired threaded hole on the skull clamp. Turn knob clockwise to tighten. Pressure applied to the Swivel Adapter in the locked position should not yield any movement in any direction.

Figure 2

Figure 3 Lock clamp in the desired position by depressing on the locking lever to its locked position. For supine and prone positions, skull clamp is positioned perpendicular to the floor and parallel for procedures requiring the sitting position. **Figure 4**

Assembly Functional Check

The following components on the M-1525 adjustable base unit must be checked before each use to ensure the proper unit functioning during each procedure.

- Pull Locking Lever to release.
- Remove the Hex Wrench from its end bracket location by turning it counterclockwise.
- With the aid of the Hex Wrench, unlock the Tension Adjustment Knob by turning the socket screw located at the center of the Knob counterclockwise.
- Turn the Tension Adjustment Knob clockwise to tighten tension and counterclockwise to loosen tension on the Locking Lever.
- Test tension on the Locking Lever by depressing and releasing Lever. Proper Lever tension is indicated by an effortless, yet firm lever movement.
- Test the overall stability of the unit, insert it into the operating room table and depress Locking Lever. Apply pressure on the Transitional Member ensuring no movement is observed.
- If movement is detected, the unit must be serviced immediately.
- Lock the Tension Adjustment Knob by turning the Hex Wrench clockwise in the socket screw.
- Return the Hex Wrench to its receptacle.

The following components on the M-1500 skull clamp must be checked before each use to ensure the proper unit functioning during each procedure.

The internal mechanism on the locking ratchet must be tested by first removing the ratchet arm from the skull clamp base. Apply upward pressure on the locking ratchet while pointing it to the floor to identify any movement on this component. If the release mechanism moves 1/16" or more, the unit must be serviced.

To ensure proper locking of the rocker arm, rotate it various times, lock and unlock by turning the index knob clockwise and validate the engagement of the locking mechanism on the rocker arm. If the rocker arm does not lock as directed, unlock the index knob, turn the rocker arm slightly and the locking mechanism should engage as expected. If locking mechanism fails, send unit for service.

The force gauge should be inspected by applying firm pressure on the skull pin receptacle until flush to the threaded area, hold and then release until pressure indicator returns to its resting position. If the pressure indicator is not flush to the pressure gauge knob, unit must be sent in for service. **Figure 5**

Transitional Member Removal

- Depress Locking Lever.
- Turn Retaining Screw counterclockwise until the Transitional Member can be easily removed.
- Insert replacement Transitional Member into its receptacle.
- Finger-tighten Retaining Screw until Transitional Member is securely in place and cannot slide out accidentally.

Skull Pin Insertion

- Push skull pins into their respective skull clamp locations until fully inserted. Skull pins must be fully seated so that there are no gaps between the shoulder of the skull pin and the receptacle into which it has been inserted.
- Remove protective skull pin sleeve before use. **Figure 6**

Patient Fixation

- Pull locking ratchet and adjust ratchet arm to desired width.
- Position the patient's head inside skull clamp ensuring that pins on the rocker arm are equidistant from the midline of the head.
- Fix rocker arm in place by rotating the index knob clockwise until arrows are aligned on skull clamp base and knob.
- Push ratchet arm toward the skull ensuring skull pins enter the patient's scalp securely at a perpendicular angle.
- Adjust fixation force by turning Force Gauge clockwise to desired strength in pounds (See Warnings Section) indicated by the torque screw on the gauge.
- Turn the Force Gauge to desired fixation strength, then back off one half turn.
- **Failure to properly position the patient and to fully tighten all adjustment portions of this device may result in patient injury.**
- Prior to pin impingement, turn the torque knob until the engagement line of the piston is visible, just beyond the clamp base. Failure to do this may cause difficulty in disengaging the skull pins after completion of the procedure.

Warranty

CareFusion offers a lifetime guarantee on every surgical device bearing the V. Mueller brand name (unless otherwise noted) to be free of functional defects in workmanship and materials when used normally for its intended surgical purpose. Any V. Mueller device proving to be defective will be replaced or repaired at no charge.

Repair Service

Regardless of age, if any V. Mueller device needs service, return it to an authorized repair service center. For repairs outside the U.S., please contact your local distributor.

Note: All devices being returned for maintenance, repair, etc. must be cleaned and sterilized per these instructions for use prior to shipment.

Contact Information

CareFusion

1500 Waukegan Rd
McGaw Park, IL 60085 U.S.A.
800-323-9088
www.carefusion.com

For domestic inquiries email:
GMB-VMueller-Cust-Support@carefusion.com

For international inquiries email:
GMB-SIT-International-Team@carefusion.com

Other Resources: To learn more about sterilization practices and what is required of manufacturers and end users, visit www.aami.org, www.aorn.org or www.iso.org.

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